

AdvanCE®
capsule endoscope
delivery device



listening... and delivering solutions®

Direct placement of video capsules...

The *AdvanCE[®] capsule delivery system* allows for direct endoscopic placement of video capsules* in patients who are unable to swallow or pass the capsule through the pylorus.

The AdvanCE[®] system...

- Facilitates endoscopic delivery of video capsules into the stomach or duodenum
- Minimizes wasted video capsules due to lost battery life spent in the esophagus or stomach
- Enables small intestinal surveillance in challenging patients with dysphagia, pill phobia, gastroparesis, anatomical abnormalities, and more

For more information on the AdvanCE[®] delivery system, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	units/box
00711144	2.5	180	3

Distributed by:



*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

*For use with the PillCam SB video capsule (not included).

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

760398 Rev. A