

# Balloon Enteroscopy



## ERCP-Cannula





Product Code	Order No.	Characteristics	Length	Accepted Guide Wire	∅ Working Channel	Unit
CAN1-B3-18-260-35	501856	Tapered, radiopaque marking, 1 lumen	260 cm	.035"	2,0 mm	2

## Papillotomes

Handle with OL power connection



Product Code	Order No.	Characteristics	Cutting Length	Accepted Guide Wire	∅ Working Channel	Unit
<b>AXS_tome</b> PAP1-D2-30-21-OL	501853	Micro-tip, 2 lumen,  , anatomical 3D orientation	30 mm	.021"	2,8 mm	2
PAP1-N2-20-21-OL	501902	B II, 2 lumen, 	20 mm	.021"	2,8 mm	2

# Balloon Enteroscopy

## EndoBite™ \_ Biopsy forceps



Product Code	Order No.	Characteristics	∅ Jaws	Length	∅ Working Channel	Unit
--------------	-----------	-----------------	--------	--------	-------------------	------





BIO1-C4-18-260	502006	Oval jaws, mono pull-wire, coated	1,8 mm	260 cm	2,0 mm	10
----------------	--------	-----------------------------------	--------	--------	--------	----

## PolySnare™ \_ Polypectomy snare



Product Code	Order No.	Characteristics	∅ Snare	Length	Handle	∅ Working Channel	Unit
--------------	-----------	-----------------	---------	--------	--------	-------------------	------



POL1-B1-30-18-260-OL	501208	Convex, 	30 mm	260 cm	OL power connection	2,0 mm	10
----------------------	--------	---	-------	--------	---------------------	--------	----


# Balloon Enteroscopy

## Injection needle



Product Code	Order No.	Characteristics	Needle ø   Length	Length	ø Working Channel	Unit
INJ1-A1-07-5-18-260	500753	Handle locking mechanism, distal metal hub, suitable for Histoacryl application	0,7 mm   5 mm	260 cm	2,0 mm	10

## Stone extraction baskets



Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
BAS1-C1-20-18-260	500025	Half twisted	20 mm	260 cm	2,0 mm	2
BAS1-C2-25-23-260	500026	Half twisted	25 mm	260 cm	2,8 mm	2



Marked products are not suitable for emergency lithotripsy.


Errors excepted.

# Balloon Enteroscopy

## Cleaning brush



Product Code	Order No.	Characteristics	ø Bristle	Length	ø Working Channel	Unit
CLE1-C2-30-16-280	500463	2 Bristles	3,0 mm	280 cm	2,0 mm	10



## Guide wire



Product Code	Order No.	Characteristics	Diameter	Length	Unit
WIR1-F1-21-500	501911	Straight, Nitinol, coated, radiopaque tip	.021"	500 cm	2

