

**medwork®**  
let's work together

Our national medical product advisors  
and international partners are pleased  
to give you detailed information about our  
ERCP system solution **easyCure®**

Feel free to contact us!



medwork medical products and services GmbH  
Dr.-Schätzel-Straße 30 | 91315 Höchstädt/Aisch | Germany  
Tel.: +49 (0) 91 93 / 6 39 19-0 | Fax: +49 (0) 91 93 / 6 39 19-99  
contact@medwork.eu | www.medwork.eu

# EndoTherapy

**medwork®**

Innovation  
Fascination  
Precision



**easyCure®**  
The Therapeutic ERCP System Solution

# easycure®

The Therapeutic ERCP System Solution

**easycure** is the innovative medwork solution for therapeutic ERCP. The special design of **easycure** instrumentation now enables you to organise ERCP processes very easily and reliably.

## easy

The **easycure** system sets new standards with its highly simplified change of instruments. You cannot only load the **easycure** instruments easily and reliably but also remove them from the securely fixed guide wire quickly by hand.

## secure


**easycure** is the secure, gentle ERCP technique. The guide wire is prevented from slipping out of the biliary system by continuous fixation during the intervention. In this way you benefit from much improved guide wire control. Even complex stone extractions and prosthesis implantations can be performed reliably and efficiently with the **easycure** system.

## easycure

The **easycure** instrumentation meets the highest standards of hygiene and functionality for ERCP. All the components in the **easycure** system are optimally matched to one another and allow you to ensure maximum efficiency in the workflow. In addition, excellent functionality and compatibility save you complicated supplementary devices and unnecessary capital investments.


### **easycure** Guide Wire Baskets in the BAS Series

With guide wire stiffening the lithotripsy-compatible **easycure** guide wire baskets can be positioned much more efficiently and are provided with a distal guide wire port. On account of the unique integrated arrangement of the guide wire lumen the papilla can be probed very gently.




### **easycure** Stone Extraction Balloons in the BAL Series

The **easycure** stone extraction balloons are also provided with a distal guide wire port. In addition, you can administer contrast media very smoothly and precisely through a separate third lumen.




### **easycure** ERCP-Cannulas in the CAN Series

The two guide wire ports of the **easycure** ERCP-cannulas enable you to flexibly probe the papilla with and without guide wire stiffening.



### **easycure** Papillotomes in the PAP Series

For customised papillotomy requirements you have various versions of our **easycure** papillotomes at your disposal. In addition to monofilament and braided cutting wires in lengths of 20 mm and 30 mm you can also take advantage of our needle papillotome. Two guide wire ports considerably facilitate the introduction and removal of instruments.



### **easycure** Guide Wires in the WIR Series

Owing to their length of 200 cm **easycure** guide wires make it much easier for you to perform ERCP processes and substantially reduce the infection risk. Furthermore, by using the 0.021" and 0.035" diameter **easycure** guide wires you can exploit the entire instrumental flexibility for ERCP to full advantage.

Benefit from the innovative **easycure** solutions for your interventional ERCP requirements!

### **easycure** Endoprosthesis Introducer Sets in the PRO Series


One or more endoprostheses can be placed very easily and reliably with the aid of **easycure** endoprosthesis introducer sets. Moreover, the radiopaque tip of the guide catheters make it easier for you to position endoprostheses with precision. You can take advantage of prostheses with a very wide range of lengths in the diameters 8.5 F and 10 F.



With its flexible properties and simple use the **easycure** adaptor makes it easier for you to integrate it perfectly into your customised ERCP processes.

### **easycure** Adaptor - SPE1-Y7

Our **easycure** adaptor provides fast and easy fixation of up to two guide wires in a stable position during ERCP. All **easycure** instruments can be introduced smoothly and moved freely inside the endoscope.



Intracatheter  
Fascination  
Precision